



# Technical Data Sheet



## nitrylex® classic (blue)

### PRODUCT DESCRIPTION

|                   |   |
|-------------------|---|
| Type of the glove | Non-sterile, powder free, examination and protective glove for single use |
| Material          | Nitrile   |
| Donning powder    | -   |
| Colour            | Blue  |
| Shape             | Ambidextrous, gloves fitting either hand                                  |
| Cuff              | Beaded  |
| External surface  | Microtextured + fingertip textured, polymerized                           |
| Internal surface  | Polymerized + chlorinated   |
| Packaging         | 10 x 100 pcs 10 x 200 pcs   |

### PRODUCT REFERENCES

| SIZE / REFERENCE NUMBER | XS | RD30019001 | XS | RD30096001 |
|-------------------------|----|------------|----|------------|
|                         | S  | RD30019002 | S  | RD30096002 |
|                         | M  | RD30019003 | M  | RD30096003 |
|                         | L  | RD30019004 | L  | RD30096004 |
|                         | XL | RD30019005 | XL | RD30096005 |

### Dimensions

#### Length [mm]

Minimum

#### Width [mm]

#### Thickness

(single wall)

[mm]

Minimum

### Elongation at

break [%]

Minimum

### Force at

break[N]

Minimum

| Size          | XS<br>(5-6) | S<br>(6-7) | M<br>(7-8) | L<br>(8-9) | XL<br>(9-10) |
|---------------|-------------|------------|------------|------------|--------------|
|               | 240         | 240        | 240        | 240        | 240          |
|               | ≤80         | 80<br>±10  | 95<br>±10  | 110<br>±10 | ≥110         |
| Middle finger |             |            | 0,06       |            |              |
| Palm          |             |            | 0,05       |            |              |
| Cuff          |             |            | 0,04       |            |              |
| Before ageing |             |            | 500        |            |              |
| After ageing  |             |            | 400        |            |              |
| Before ageing |             |            | 6,0        |            |              |
| After ageing  |             |            | 6,0        |            |              |

### MANUFACTURING AND SAFETY STANDARDS

|                                |   |   |
|--------------------------------|---|---|
| AQL                            | Manufacturing final release: G-I inspection level AQL 1.0 in accordance with ISO 2859-1   |   |
| Protein content                | N/A   |   |
| CE classification              | Class I – Medical Device (Council Directive 93/42/EEC)  | Category III – Personal Protective Equipment (Regulation (EU) 2016/425)         |
| Compliances                    | EN 455-1, EN 455-2, EN 455-3, EN 455-4<br>EN ISO 15223-1<br>EN 1041<br>EN ISO 13485   | EN ISO 374-1 (Type B), EN 374-2, EN 374-4, EN ISO 374-5<br>EN 16523-1<br>EN 420 |
| Viral test                     | Test in accordance with ASTM F1671 & ISO 16604  |   |
| Cytostatics permeation         | Test in accordance with ASTM D6978  |   |
| Chemical substances permeation | Test in accordance with EN 16523-1  |   |
| Food contact                   | Declaration of conformity for food contact in accordance with Regulation (EC) No 1935/2004 and with Commission Regulation (EU) No 10/2011 and Overall Migration Test in accordance with Commission Regulation (EU) No 10/2011 |   |
| Shelf life                     | 3 years   |   |

### STORAGE

|                     |   |  |
|---------------------|---|--|
| Storage instruction | Keep out of direct sunlight. Store in a cool, dry place in temperature 5-35° C. Keep away from sources of ozone and ignition.<br>Issue date: 30.05.2017<br>Update: 19.11.2019 | Prepared by: Wojciech Hercka – Product Documentation Manager<br>Approved by: Katarzyna Nowakowska – Senior Product Manager |
|---------------------|---|--|



## nitrylex® classic

The instruction below should be used in conjunction with detailed information on the packaging.

### Short description of the product

Nitrile examination and protective gloves, powder-free, non-sterile for disposable use

### Full description of the product

|                       |  |
|-----------------------|--|
| Raw material          | : nitrile  |
| External surface      | : bisque with fingertip textured, polymerized      |
| Internal surface      | : polymerized + chlorinated                        |
| Cuff                  | : beaded   |
| Colour                | : blue/white/violet                                |
| Shape                 | : ambidextrous, fitting to the right and left hand |
| Size range            | : XS (5-6), S (6-7), M (7-8), L (8-9), XL (9-10)   |
| AQL                   | : 1.0  |
| Quantity in packaging | : 50/100/200 pcs. by weight                        |
| Shelf life            | : 3 years (from the date of manufacturing)         |


### Storage instructions

It is recommended to store the gloves in dry place, in the temperature of 5-35°C and to protect them against direct sunlight and fluorescent light. Recommended relative humidity in the room where the gloves are stored is 60 ±20%.

Keep the gloves in a distance of not less than 1m from heating devices, sources of fire and ozone.

Do not keep in direct vicinity of solvents, oils, fuels and lubricants.

### Food contact

Gloves are marked with food contact symbol  and comply with the requirements of Regulation (EU) No 10/2011, European Regulation (EC) No 1935/2004 and with Regulation (EC) No 2023/2006 on Good Manufacturing Practice. Gloves are suitable for handling any type of food and have been tested for Overall Migration Test acc. EN 1186:

| Extraction conditions<br>(tested for 2 h in 40°C) | Analysis results<br>[mg/dm <sup>2</sup> ] | Test Result<br>(limit < 10 mg/dm <sup>2</sup> ) |
|---|---|---|
| 3% acetic acid                                    | 1,1                                       | Pass  |
| 10% ethanol                                       | <1  | Pass  |
| Olive oil   | <3  | Pass  |

### MDD classification & compliance

Gloves are classified as class I Medical Device as per Annex IX of the Council Directive 93/42/EEC and comply to standards:

EN 455-1:2000, EN 455-2:2015, EN 455-3:2015, EN 455-4:2009, EN ISO 15223-1:2016, EN 1041:2008+A1:2013.

### PPE classification & compliance

Gloves are category III Personal Protective Equipment as per Annex I of the Regulation 2016/425 and comply to standards:


EN 420:2003+A1:2009, EN ISO 374-1:2016 (Type B), EN 374-2:2014, EN 16523-1:2015, EN 374-4:2013, EN ISO 374-5:2016.

Declaration of Conformity can be found under below web address:

<http://mercatormedical.eu/produkty/rekawice/diagnostyczne/nitrylex-classic>

**Notified Body 2777**  
**responsible for EU Type**  
**Examination (Module B)**  
**and Module C2 On-going**  
**Conformity:**

**Satra Technology Europe Ltd**  
Bracetown Business Park, Clonee  
Dublin 15, Dublin, Ireland



### Intended use

These are non-sterile examination and protective gloves for single use, intended for use in medical field to: protect patient and user from cross-contamination, conducting medical examinations, diagnostic and therapeutic procedures and for handling medical contaminated material. Gloves are classified as Medical Devices Class I and as a Personal Protective Equipment category III. Their design and labelling corresponds to the requirements of the European Medical Device Directive 93/42/EEC and the European Regulation 2016/425 on Personal Protective Equipment. Gloves should be used solely according to their intended application.

### Precautions and indications for use

Dry hands before putting the gloves on. Before usage, inspect the gloves for any defect or imperfections. Use at least 1 pair of gloves for one patient and one procedure, these are disposable gloves. Do not let chemical substances get under the gloves through the cuff. If a chemical substance reaches the skin, wash it away immediately with plenty of water with soap. If the gloves get punctured, torn or broken during their use, take them off and put on the new ones. Avoid using gloves dirty in the inside as they may cause irritation leading to skin inflammation or more serious damages. The gloves should not be used in contact with open fire and to protect against any sharp tools. The gloves are not intended for welding, electric shock protection, ionizing radiation or from the effect of hot or cold objects.

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.

The chemical penetration resistance has been assessed under laboratory conditions from samples taken from the palm only (except in case where glove is equal to or over 400 mm – where the cuff is tested also) and relates only to the chemical tested and to the tested specimen. It can be different if the chemical is used in a mixture.

It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on the temperature, abrasion and degradation.

When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.

Gloves are suitable for special purposes as they are examination gloves where risk of injury to the wrist is considered to be minimal, gloves are shorter than EN 420 min. length requirement.

### Components / hazardous components

Some gloves may contain components known to be a possible cause of allergy for person allergic to them, who may develop contact irritation and/or allergic reaction. In case of an allergic reaction, seek medical assistance immediately.

### Disposal

Used gloves can be contaminated with contagious or other hazardous substances. They should be disposed of in accordance with local regulation. Gloves should be buried or burned under controlled conditions.

### Manufacturer

MERCATOR MEDICAL S.A.

ul. H. Modrzejewskiej 30

31-327 Cracow, Poland

[www.mercatormedical.eu](http://www.mercatormedical.eu)





### Permeation performance levels as per EN ISO 374-1:2016

• Level 1 > 10 min • Level 2 > 30 min • Level 3 > 60 min • Level 4 > 120 min • Level 5 > 240 min • Level 6 > 480 min

| Test results acc. to EN 16523-1:2015 |       | EN 374-4:2013   | Test results acc. to EN 16523-1:2015 |       | EN 374-4:2013   |
|--------------------------------------|-------|-----------------|--------------------------------------|-------|-----------------|
| Chemical                             | Level | Degradation [%] | Chemical                             | Level | Degradation [%] |
| *4% Chlorhexidine Digluconate        | 6     | 19.0            | 30% Hydrogen Peroxide (P)            | 2     | 22.8            |
| 40% Sodium Hydroxide (K)             | 6     | -42.9           | 1.5% Methanol in water               | 6     | 21.9            |
| 10-13% Sodium Hypochlorite           | 6     | 14.7            | 25% Ammonium Hydroxide (O)           | 1     | -52.0           |
| 50% Sulphuric Acid                   | 6     | -20.5           | 3% Povidone-iodine                   | 6     | 33.7            |
| 10% Acetic Acid                      | 4     | 66.7            | 10% Sodium Percarbonate              | 6     | 15.4            |
| 5% Ethidium Bromide                  | 6     | 3.4             | 50% Glutaraldehyde                   | 6     | 27.4            |
| 37% Formaldehyde (T)                 | 3     | 5.0             | 0.1% Phenol                          | 6     | 33.8            |

\*Permeation rate 7µg/cm²/min, EN 374-4:2013 Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical.

| Test acc. To EN 374-2:2014 – Level 2 (ISO 2859) |        | Test acc. To EN ISO 374-5:2016      |      |
|---|--------|-------------------------------------|------|
| Performance level                               | AQL    | Protection against bacteria & fungi | Pass |
| Level 3   | < 0.65 | Protection against viruses          | Pass |
| Level 2   | < 1.5  |                                     |      |
| Level 1   | < 4.0  |                                     |      |



### Symbols used on the packaging

|  |  |  |  |  |   |
|--|--|--|--|--|---|
|  | Do not re-use / gloves are intended for single use   |  | Non-sterile gloves   |  | Powdered gloves   |
|  | Do not use, if package is damaged  |  | Keep away from solar and fluorescent light                                     |  | Powder free gloves  |
|  | Keep away from moisture, store in a dry place  |  | Temperature limitation / gloves store in temperature 5-35°C                    |  | Presence of polymer coating on the inner surface of the glove |
|  | Raw material – natural rubber latex  |  | Keep away from ozone   |  | Presence of cosmetic coating on the inner surface of glove    |
|  | Catalogue number   |  | Lot / batch number   |  | Gloves with incorporated singlet oxygen layer.                |
|  | EU Authorized Representative, symbol should be accompanied by name and address of Authorized Representative                  |  | Expiry date  |  | Presence of external texture on the glove                     |
|  | Marking of gloves protecting against bacteria and fungi.   |  | Gloves protecting against chemical dangers with digit literal odes             |  | Gloves made from nitrile                                      |
|  | Marking of gloves protecting against viruses, bacteria and fungi.  |  | Antistatic gloves  |  | Gloves made from vinyl  |
|  | Marking o type A chemical resistant gloves. Six tested chemicals shall be identified by their code letter under pictogram.   |  | Date of manufacture  |  | Gloves made from neoprene                                     |
|  | Marking o type B chemical resistant gloves. Three tested chemicals shall be identified by their code letter under pictogram. |  | Manufacturer, symbol should be accompanied by name and address of Manufacturer |  | Gloves made from polyisoprene                                 |
|  | Marking o type C chemical resistant gloves. One tested chemicals shall be identified by their code letter under pictogram.   |  | Consult instructions for use   |  | 50 gloves by weight   |
|  | Protective glove against mechanical risk (if applicable accompanied by 4 digit code of relevant performance levels)          |  | Package made from paper, qualify for recycling                                 |  | 100 gloves by weight  |
|  | Food contact symbol (article is suitable for food contact, for details check the instruction for use)                        |  | Package is treated as municipal waste  |  | 200 gloves by weight  |
|  | Indicates compliance with the requirements of Russian market   |  | Indicates compliance with the requirements of Ukrainian market                 |  | Additional information on inner side of package               |










## ■ HOW TO PUT THE GLOVES ON?



## ■ HOW TO TAKE THE GLOVES OFF?





**EU DECLARATION OF CONFORMITY**

Manufacturer: **MERCATOR MEDICAL S.A.**  
UL. H.MODRZEJEWSKIEJ 30  
31-327 KRAKÓW, POLSKA

Declares under its sole responsibility that non-sterile examination and protective gloves:

| Brand                    | Type   | Sizes                | Reference Numbers   |
|--------------------------|--|----------------------|---|
| <b>nitrilex® classic</b> | nitrile, powder-free, blue, for single use   | XS (5-6) - XL (9-10) | a'100: RD30019001-05<br>a'200: RD30096001-05                        |
|                          | nitrile, powder-free, white, for single use  | XS (5-6) - XL (9-10) | a'50: RD30174001-05<br>a'100: RD30143001-05<br>a'200: RD30097001-05 |
|                          | nitrile, powder-free, violet, for single use | XS (5-6) - XL (9-10) | a'100: RD30169001-05<br>a'200: RD30168001-05                        |

classified as medical device class I according to Annex IX of the Council Directive 93/42/EEC meet the essential requirements of Annex I of the Council Directive 93/42/EEC amended by the Directive 2007/47/EC and comply with the European harmonized standards: EN 455, EN ISO 15223-1, EN 1041. Conformity assessment procedure performed according Annex I and Annex VII of the Council Directive 93/42/EEC amended by the Directive 2007/47/EC.

The products described above are also classified as Personal Protective Equipment Category III and comply with Regulation (EU) 2016/425 of the European Parliament and the Council of 9 March 2016 on Personal Protective Equipment and resolution of the Council Directive 89/686/EEC and European standards: EN 420:2003+A1:2009, EN ISO 374-1:2016, EN 374-2:2014, EN 16523-1:2015, EN 374-4:2013, EN ISO 374-5:2016.

The products described above are identical to the Personal Protective Equipment, which is the subject to the EU Type Examination (Module B) under certificate No. 2777/10015-03/E17-01 issued by notified body:

**Satra Technology Europe Limited (2777)**

Bracetown Business Park, Clonee, Dublin 15, Dublin, Ireland

and are subject to the conformity to type procedure based on the internal production control plus supervised product checks at random intervals (Module C2) under surveillance of the notified body:

**Satra Technology Europe Limited (2777)**

Bracetown Business Park, Clonee, Dublin 15, Dublin, Ireland

Date and place of issue:  
30.09.2019, Kraków

**MERCATOR MEDICAL S.A.**  
ul. Heleny Modrzejewskiej 30, 31-327 Kraków  
tel. 12 66 55 400, fax 12 66 55 415  
Rejestracja: Sąd Rejonowy dla Krakowa - Śródmieścia w Krakowie,  
XI Wydział Gospodarczy KRS, KRS: 0000036244  
Kapitał zakładowy (w całości wpłacony): 10.589.100 PLN  
NIP: 677-10-36-424, REGON: 350967107  
Numer BDO: 000056063  
-6-

Signed on the behalf of the Manufacturer:



Wojciech Hercka  
Product Documentation Manager

# EU Type-Examination Certificate

## Certificate number: 2777/10015-03/E17-01

This EU Type-Examination Certificate covers the following product group(s) supported by testing to the relevant standards/technical specifications and examination of the technical file documentation:  
Following the EU Type-Examination this product group has been shown to satisfy the applicable essential health and safety requirements of Annex II of the PPE Regulation (EU) 2016/425 as a Category III product.

### Product Reference

Nitrylex classic  
Nitrylex one by one

### Description

Disposable powder free nitrile glove available in white, blue and violet blue. Non-Sterile

Sizes: 5/6 (XS) to 9/10 (XL)

### Classification:

#### EN ISO 374-1:2016 (Type B)

|                               | Level | EN374-4:2013 Degradation % |
|-------------------------------|-------|----------------------------|
| *4% Chlorhexidine Digluconate | 6     | 19.0                       |
| 40% Sodium Hydroxide (K)      | 6     | -42.9                      |
| 10-13% Sodium Hypochlorite    | 6     | 14.7                       |
| 50% Sulphuric Acid            | 6     | -20.5                      |
| 10% Acetic Acid               | 4     | 66.7                       |
| 5% Ethidium Bromide           | 6     | 3.4                        |
| 37% Formaldehyde (T)          | 3     | 5.0                        |
| 50% Glutaraldehyde            | 6     | 27.4                       |
| 0.1% Phenol                   | 6     | 33.8                       |
| 30% Hydrogen Peroxide (P)     | 2     | 22.8                       |
| 1.5% Methanol in water        | 6     | 21.9                       |
| 25% Ammonium Hydroxide (O)    | 1     | -52.0                      |
| 3% Povidone-iodine            | 6     | 33.7                       |
| 10% Sodium Percarbonate       | 6     | 15.4                       |
| * Permeation rate 7µg/cm²/min |       |                            |

#### Level

#### EN374-4:2013 Degradation %

#### EN ISO 374-5: 2016

Protection against Bacteria & Fungi  
Protection against viruses

Pass  
Pass

### Standards/Technical specifications applied:

EN 420: 2003+A1: 2009; EN ISO 374-1:2016; EN ISO 374-5:2016


### Technical reports/Approval documents:

SATRA: SPC0216113/1327/SMcD/RS, SPC0216113/1327, PRC0250570/1640/SPT, CHM0248297/1630/EN/E, CHM0248297/1630/EN/D/Issue 2, CHM0257198/1719/SMcD/A, CHM0257198/1719/SMcD/B, CHM0257198/1719/SMcD/C, CHM0257198/1719/SMcD/D

Signed on behalf of SATRA:



Hannah Coe



Jacquie Glasspool

Date of issue: 18/09/2019

Expiry date: 21/04/2023

# TERMS AND CONDITIONS

The following conditions apply in addition to SATRA's standard terms and conditions of business and those given in the current certification agreement.

The certificate holder is licensed to mark the products detailed within this certificate in accordance with Annex V (Module B) of the Regulation (EU) 2016/425 of the European Parliament and of the council of 9th March 2016 on personal protective equipment once you have drawn up an EU declaration of product conformity.

Please note:

1. Where the product is classified as category III then CE Marking of production is reliant on current compliance with Regulation 2016/425 module C2 or Module D. (Except that specifically produced to fit an individual user).
2. Full details of the scope of the certification and product(s) certified are contained within the manufacturer's technical documentation.
3. Where a translation of this certificate exists, the English language version shall be considered as the authoritative text.
4. Certification is limited to production undertaken at the sites listed in the manufacturers technical documentation.
5. Ongoing manufactured product shall be consistent with the product(s) certified and listed on this certificate.
6. The Manufacturer shall inform SATRA of any changes to the certified product or technical documentation.
7. Where results obtained during type testing are within the budget of uncertainty when compared to the pass requirement, classification or performance level, then it is the responsibility of the manufacturer to ensure that the factory production control and manufacturing tolerances are such that the product placed on the market meets with the stated requirements, classifications or performance levels.
8. This certificate shall be kept together with the relevant technical documentation in a safe place by the client named on this certificate. Production of this certificate and other documentation may be required by a representative of the EC member state government.
9. This certificate relates only to the condition of the testable items at the time of the certification procedure and is subject to the expiry date shown.
10. SATRA reserves the right to withdraw this certificate if it is found that a condition of manufacture, design, materials or packaging have been changed and therefore no longer comply with the requirements of Regulation 2016/425.